



PC1700 2004 / 0 0 2 0 2 1



INVESTOR IN PEOPLE

The Patent Office
Concept House
Cardiff Road
Newport
South Wales
NP10 800

**PRIORITY
DOCUMENT**
SUBMITTED OR TRANSMITTED IN
COMPLIANCE WITH RULE 17.1(a) OR (b)

REC'D 01 JUL 2004

WIPO

PC1

I, the undersigned, being an officer duly authorised in accordance with Section 74(1) and (4) of the Deregulation & Contracting Out Act 1994, to sign and issue certificates on behalf of the Comptroller-General, hereby certify that annexed hereto is a true copy of the documents as originally filed in connection with the patent application identified therein.

I also certify that the attached copy of the request for grant of a Patent (Form 1/77) bears an amendment, effected by this office, following a request by the applicant and agreed to by the Comptroller-General.

In accordance with the Patents (Companies Re-registration) Rules 1982, if a company named in this certificate and any accompanying documents has re-registered under the Companies Act 1980 with the same name as that with which it was registered immediately before re-registration save for the substitution as, or inclusion as, the last part of the name of the words "public limited company" or their equivalents in Welsh, references to the name of the company in this certificate and any accompanying documents shall be treated as references to the name with which it is so re-registered.

In accordance with the rules, the words "public limited company" may be replaced by p.l.c., plc, P.L.C. or PLC.

Re-registration under the Companies Act does not constitute a new legal entity but merely subjects the company to certain additional company law rules.

Signed

Andrew Gurney

Dated 24 June 2004

Patents Form 1/77
Patent Act 1977
(Rule 16)

THE PATENT OFFICE
K
16 JUN 2003
NEWPORT



16 JUN 03 E815176-1 B11618
P01/7700 0.00-0313828.6

Request for grant of a patent
(See the notes on the back of this form. You can also get an explanatory leaflet from the Patent Office to help you fill in this form)

The Patent Office
Cardiff Road
Newport
South Wales
NP10 8QQ

1. Your reference

B05001

2. Patent application number
(The Patent Office will fill in this part)

0313828.6

16 JUN 2003

3. Full name, address and postcode of the or of each applicant (underline all surnames)

BRIAN DAVID OWEN OWEN-SMITH
48, WESTGATE
CHICHESTER
WEST SUSSEX
PO19 3EU

8653883001

Patents ADP number (if you know it)
If the applicant is a corporate body, give the country/state of its incorporation

4. Title of the invention

SALIVARY URATE FOR DIAGNOSIS OF PRE-ECLAMPSIA

5. Name of your agent (if you have one)

"Address for service" in the United Kingdom to which all correspondence should be sent (including the postcode)

KILBURN + STRODE
20 REGENT ST
LONDON
WC1R 4PJ

Patents ADP number (if you know it)

6. If you are declaring priority from one or more earlier patent applications, give the country and the date of filing of the or of each of these earlier applications and (if you know it) the or each application number

Country

Priority application number
(if you know it)

Date of filing
(day / month / year)

125001

FS1A7

7. If this application is divided or otherwise derived from an earlier UK application, give the number and the filing date of the earlier application

Number of earlier application

Date of filing
(day / month / year)

8. Is a statement of inventorship and of right to grant of a patent required in support of this request? (Answer 'Yes' if:
a) any applicant named in part 3 is not an inventor, or
b) there is an inventor who is not named as an applicant, or
c) any named applicant is a corporate body. See note (d))

No

Patents Form 1/77

9. Enter the number of sheets for any of the following items you are filing with this form. Do not count copies of the same document

Continuation sheets of this form

NONE

Description

~~SIX~~ FIVE

Claim(s)

ONE

Abstract

Drawing(s)

SIX ONLY

10. If you are also filing any of the following, state how many against each item.

Priority documents

Translations of priority documents

Statement of inventorship and right to grant of a patent (Patents Form 7/77)

Request for preliminary examination and search (Patents Form 9/77)

Request for substantive examination (Patents Form 10/77)

Any other documents (please specify)

11.

I/We request the grant of a patent on the basis of this application.

Signature

B.D. Owen-Smith

Date

14th June 2003

12. Name and daytime telephone number of person to contact in the United Kingdom

01243 786688

Warning

After an application for a patent has been filed, the Comptroller of the Patent Office will consider whether publication or communication of the invention should be prohibited or restricted under Section 22 of the Patents Act 1977. You will be informed if it is necessary to prohibit or restrict your invention in this way. Furthermore, if you live in the United Kingdom, Section 23 of the Patents Act 1977 stops you from applying for a patent abroad without first getting written permission from the Patent Office unless an application has been filed at least 6 weeks beforehand in the United Kingdom for a patent for the same invention and either no direction prohibiting publication or communication has been given, or any such direction has been revoked.

Notes

- If you need help to fill in this form or you have any questions, please contact the Patent Office on 08459 500505.
- Write your answers in capital letters using black ink or you may type them.
- If there is not enough space for all the relevant details on any part of this form, please continue on a separate sheet of paper and write "see continuation sheet" in the relevant part(s). Any continuation sheet should be attached to this form.
- If you have answered 'Yes' Patents Form 7/77 will need to be filed.
- Once you have filled in the form you must remember to sign and date it.
- For details of the fee and ways to pay please contact the Patent Office.

Salivary urate for diagnosis of pre-eclampsia

DUPLICATE

Dr Brian D Owen Smith

Description

[0001] This invention relates to detection of pre-eclampsia. In particular, it relates to diagnosis of pre-eclampsia by detecting elevated levels of urate, the sodium salt of uric acid, 2,6,8, trioxypurine, in a biological fluid, saliva.

[0002] The invention is a method of using the enhanced sensitivity inherent in urate concentration in saliva, as compared with that of blood, and depends on urate as a metabolic compound with unique and predictable behaviour.

[0003] Pre-eclampsia occurs in approximately 2-4% of pregnant women and is a common cause of maternal mortality and foetal and maternal morbidity. In severe cases, early delivery may be necessary and there is therefore the risk of the child being handicapped.

[0004] Pre-eclampsia, otherwise known as 'gestational proteinuric hypertension (GPH)' is a multisystem disease of pregnancy of unknown cause. The maternal syndrome is characterised by various abnormalities: increased blood pressure, oedema, proteinuria and abnormal clotting, liver and renal function. This may lead to heart and kidney failure, intracranial haemorrhage and fits. The cause of pre-eclampsia is believed to lie in the placenta and there is evidence for a circulating endothelial cell "toxic" factor.

[0005] Oxidative stress is believed central to development of pre-eclampsia and associated with poor placental perfusion. Raised blood urate, a marker of oxidative stress, is associated with severity of pre-eclampsia and foetal outcome (Many et al, 1996). Blood urate is also a predictor of pre-eclampsia (Chappell et al 2002).

[0006] The symptoms of pre-eclampsia are generally detectable from around 28 weeks up to full term and are not usually apparent before 24 weeks. The conventional tests, as disclosed in the prior art EP 0927355B1 are for kidney failure by measuring urea in the blood or protein in the urine, and for increased maternal blood pressure.

[0007] The prior art however fails in several important aspects, blood tests are time consuming, need an experienced operator, may cause physical and psychological discomfort and could lead to dangerous needle stick injuries such as infectious hepatitis. Blood samples also require specialist storage techniques to arrest decomposition.

[0008] What is needed therefore is a method for identification of those pregnant women at risk of pre-eclampsia with a biochemical test to detect oxidative stress, that is of high sensitivity, is inexpensive, non-invasive and is both simple and quick to perform routinely or in emergency and is safer, easier, quicker, and cheaper to obtain than blood.

[0009] A test for urate is ideal since it is a powerful antioxidant and marker of oxidative stress. Further, biochemical change is likely to occur before irreversible tissue damage. But blood urate concentration is a poor discriminator of pre-eclampsia as it is also increases in normal pregnancy.

[0010] However, blood urate is a product of endogenous urate production, diet, renal excretion and gastro-intestinal excretion. In the intestine, urate undergoes bacterial degradation (enteral uricolysis) (Sorensen, 1960) and is subject to diurnal variability (Owen-Smith et al, 1981, Figure 1). At night, urine excretion of uric acid is reduced whereas intestinal elimination is increased. It has been shown that blood urate concentrations also demonstrate diurnal variation (Kanabrocki et al, 2000, Figure 2). The diurnal variation is thought primarily due to sleep (Owen-Smith, 1983, Figure 3). The net result of these changes is relative stability of blood urate so that with overproduction of cellular urate, increase in blood concentration is delayed.

[0011] Urate is a small molecule (Molecular weight 168) and widely diffused in extracellular fluid, including saliva. Salivary urate is not reabsorbed (as in the kidney) but is swallowed and degraded by bacteria in the intestine. Thus, rapid changes of surges of urate production are detected in saliva without significant change in blood urate concentration.

[0012] It has been demonstrated that concentration of urate in saliva is a more sensitive indicator of cellular urate metabolism than is found in blood. This was established in patients suffering from gout, before and after treatment with allopurinol (that prevents urate formation) and in normal subjects undertaking exercise or providing samples day and night (to assess diurnal variability in excretion) (Owen-Smith et al, 1998).

[0013] Further, a patient with pre-eclampsia provided both blood and salivary urate samples on the days before and after delivery (Figure 4). Both sets of samples were greatly raised as compared with normal non-pregnant women and a 28 weeks pregnancy. The salivary levels in pre-eclampsia showed good correlation with blood levels. The increase in saliva concentration from 28 weeks pregnancy to delivery was 180 $\mu\text{mol/l}$.

[0014] The present invention provides in one aspect a diagnosis of pre-eclampsia which method comprises measuring urate in a biological sample for comparison with normal values.

[0015] The sample is maternal salivary urate, which has not previously been used for diagnosis of pre-eclampsia in pregnancy.

[0016] The invention may also diagnose conditions of oxidative stress from other causes and may assess activity or severity of pre-eclampsia.

[0017] The ease of collection of saliva and assay of urate fulfil the criteria of cost, safety and convenience for patient and laboratory technician.

[0018] The invention is a useful adjunct to surveillance of pregnancy and pre-eclampsia.

According to the present invention a method for testing for pre-eclampsia in pregnant human females is characterised by a method which comprises measuring urate in a biological sample.

The invention will now be described by example and with reference to the following figures.

Figure 1 Daily variation of enteral uricolysis increased at night. Graph of Disintegrations/sec/mmol $^{14}\text{CO}_2$ derived from uric acid-2- ^{14}C (1.19mmol/l) injected intravenously during (A) 5 days fasting and hyperuricaemic and (B) 5 days whilst normouricaemic on a low purine diet.

Figure 2 Showing increased sensitivity and circadian variation of salivary urate concentration as compared with blood urate using serum data from Kanabrocki et al. Salivary urate correlates with serum. Serum urate (bars and left Y axis) ranges from 380-405 $\mu\text{mol/l}$ whilst saliva (line and right Y axis) is amplified from 125- 350 $\mu\text{mol/l}$ i.e. a range of 225 $\mu\text{mol/l}$.

Figure 3 Showing fall in blood urate during 42 hours of sleep deprivation (while on a low purine diet) with increase during sleep and associated fall in urine uric acid excretion.

Figure 4 Blood and saliva data for non pregnant women, a normal pregnancy of 28 weeks and 3 days blood and saliva urate data from a patient with pre-eclampsia, who was delivered on day 2, showing elevations in salivary urate.

Figure 5 Showing a plastic capped Salivette® tube and cotton swab insert. Note insert vessel, which has a perforation in its base. From 0.5-1ml of saliva collects in the apex of the tube after centrifugation.

Figure 6 Showing mean normal diurnal variation of salivary urate in 4 non-pregnant females. The variation ranges from 70-190 $\mu\text{mol/l}$.

SAMPLE COLLECTION


[0019] Before cleaning teeth or having a meal the woman removes the cotton wool swab (Figure 5) (1), from the sterile Salivette® tube (2). A Salivette® is a plastic disposable centrifuge tube with an insert to take a cotton swab that has a perforation at its base. This is not material to the invention but we envisage other methods of saliva sampling and testing to be developed such as dry test stick or dip stick coated with suitable reagents. The woman places the cotton wool swab under the tongue to moisten it and then gently chews upon it for about 45 seconds. This activates salivation within the mouth and the swab absorbs the saliva. The swab is then returned to the insert in the Salivette®, which is appropriately labelled, to ensure safe storage of the sample. It is then kept cool in a normal domestic refrigerator at approximately 4°C until it is transported and processed by the laboratory.

[0020] Upon receipt at the laboratory and after laboratory and patient documentation is collected from the tube label, the Salivette® tube with the swab insert is centrifuged at 3000 revs for 5 minutes. The clear sample of mixed saliva (3) is then ready for assay. Assay may be undertaken by any method that provides a level of accuracy within +/- 5% and allow interpretation within a set of values from 30 $\mu\text{mol/l}$ to 1000 $\mu\text{mol/l}$. However for completeness a method using a commercially available multi-tester, the SYNCHRON CX Systems MULTI Calibrator will now be described.

[0021] The sample is its own reagent and in conjunction with the SYNCHRON CX Systems MULT Calibrator, although intended for the quantitative determination of uric acid (urate) concentration in serum, plasma or urine, saliva is perfectly compatible with this system.

[0022] Uric acid reagent is used to measure the uric acid concentration by a timed end point method. Uric acid is oxidised by uricase to produce allantoin and hydrogen peroxide. The hydrogen peroxide reacts with 4- aminoantipyrine (4-AAP) and 3,5-dichloro-2-hydroxybenzene sulfonate (DCHBS in a reaction catalysed by peroxidase to produce a coloured product.

[0023] The Synchron CX System automatically proportions the appropriate sample and reagent volumes into a cuvette. The ratio used is one part sample to 25 parts reagent. The system monitors the change in absorbance at 920 nanometers. This change in absorbance is directly proportional to the concentration of uric acid in the sample and is used by the Synchron cx system to calculate and express the uric acid concentration.



Examples

[0024] In the assessment of sensitivity of salivary urate as compared with blood (Figure 2) the results of 2 hourly urate concentrations over 24 hours were compared with the serum urate data of Kanabrocki et al. Salivary urate was correlated with the diurnal variation of serum urate, increased at night. The 24-hour range for serum was from 380 $\mu\text{mol/l}$ to 405 $\mu\text{mol/l}$ i.e. a variation of only 25 $\mu\text{mol/l}$ whilst the results for saliva were amplified from 125 $\mu\text{mol/l}$ to 350 $\mu\text{mol/l}$ i.e. a range of 225 $\mu\text{mol/l}$

[0025] Four non-pregnant women (Figure 6) gave salivary urate samples to demonstrate diurnal variation. There was increase in the early hours. This is must be taken into consideration in interpretation of test results and also changes in normal variation due to pregnancy. The average of the 8am to 6pm test results (76 $\mu\text{mol/l}$) show a significant lower value set than the pre-eclamptic (325 $\mu\text{mol/l}$) (fig 4).

Claims

- 1 A method for diagnosis of pre-eclampsia which comprises measuring urate in a biological sample.
2. A method as claimed in claim 1, wherein the sample is from maternal saliva.
3. A method as claimed in the above claims, wherein the urate is measured by means of a timed end point method.
4. A method as claimed in claim 1 or 2, wherein the urate is measured by means of a dry test method.
5. A method as claimed in any of the above claims, wherein the urate is measured by means of a dip stick method.
6. The use of urate levels in the saliva as an indicator of pre-eclampsia.
7. A diagnostic kit for performing the method of any one of the above claims, wherein the kit comprises reagents required to determine the level of the urate being measured.
8. A test equipment for performing the method in claims 1 and 2 of, wherein the test equipment comprises reagents required to determine the level of the urate in saliva being measured to determine pre-eclampsia.

REFERENCES

Many A., Hubel CA., Roberts JM. (1996) Hyperuricaemia and xanthine oxidase in preeclampsia, revisited. *American Journal of Obstetrics and Gynaecology* **174**(1), 288-291

Chappell LC., Seed PT., Briley AL., Kelly FJ., Hunt BJ., JH., Charnock-Jones S., Mallet AI., Poston L (2002). A longitudinal study of biochemical variables in women at risk of preeclampsia. *American Journal of Obstetrics and Gynaecology*, **186**, 127-136.

Sorensen LB. The elimination of uric acid in man studied by means of C¹⁴-labelled uric acid (1960) *Scandinavian Journal of Clinical Laboratory Investigation*; supplement **54**.

Owen-Smith BD, & Whyman A. (1981). Diurnal and nocturnal variations of enteral uricolysis during fasting and refeeding. (Abstract) *Annals of Rheumatic Diseases*, **40**, 523-524.

Kanabrocki EL., Third J., Ryan MD., Nemchausky BA., Shirazi P., Scheving LE., Hines E Jr., McCormick JB (2000) Circadian relationships of serum uric acid and nitric oxide. *Journal of the American Medical Association*; **283**, 2240-2241.

Owen-Smith BD (1983). Significance of variations of uric acid elimination *Annals of Rheumatic Diseases*; **42**, Supplement 87.

Owen-Smith BD., Quiney, J & Read J (1998). Salivary urate in gout, exercise and diurnal variation. *Lancet*; **351**:1932

Figure 1 Daily variation of enteral uricolysis increased at night. Graph of Disintegrations/sec/mmol $^{14}\text{CO}_2$ derived from uric acid-2- ^{14}C (1.19mmol/l) injected intravenously during (A) 5 days fasting and hyperuricaemic and (B) 5 days whilst normouricaemic on a low purine diet.




Figure 2 Comparative sensitivity of circadian variation of blood and saliva (serum data from Kanabrocki et al)




Figure 3. Effect of sleep and 42 hours sleep deprivation on blood and urine urate concentration




Figure 4. Showing blood and saliva data for non pregnant women, a normal pregnancy of 28 weeks and 3 days blood and saliva urate data from a patient with pre-eclampsia who was delivered on day 2, showing elevation in salivary urate.




Figure 5 Showing a plastic capped Salivette® tube and cotton swab insert. Note insert vessel, which has a perforation in its base. From 0.5-1ml of saliva collects in the apex of the tube after centrifugation.




Figure 6. Mean daily variation of salivary urate in 4 non-pregnant women

